



# EIA for Quantitative Determination of anti-β2-adrenergic Receptor Antibodies

## Introduction

The β-2 adrenergic receptor is a G protein-coupled receptor. β receptors have the subtypes β<sub>1</sub>, β<sub>2</sub> and β<sub>3</sub>. Autoantibodies against β<sub>2</sub> adr-r are significantly elevated in patients suffering from chronic fatigue syndrome (CFS) / myalgic encephalomyelitis (ME). β-2 adrenergic receptor (b2 adr-r) autoantibodies are present in a subset of 20-30% of chronic fatigue syndrome patients.

The CellTrend anti-β2-adrenergic-receptor-antibody-EIA is designed for the determination of antibodies (IgG) against the β2-adrenergic-receptor in serum and plasma.

## Principle of the assay

The CellTrend anti-β2-adrenergic-receptor-antibody-EIA is an antibody screening test. The β2-adrenergic-receptor has been pre-coated onto a microtiter plate. During the first incubation the anti-β2-adrenergic-receptor-antibodies of the samples are immobilised on the plate. The auto-antibodies are detected with a POD labeled anti-human IgG antibody. In the following enzymatic substrate reaction the intensity of the colour correlates with the concentration and/ or avidity of anti-β2-adrenergic-receptor-antibodies

## Performance Characteristics

### Standard curve:

5 standards between 2.5 U/ml and 40 U/ml

### cut off (Löbel et al, Brain Behav Immun 2015):

positive (95% percentile): 14.0 Units/ml

at risk (90% percentile): 8.0 – 14.0 Units/ml

### Sample materials:

Serum

### Intraassay-Precision:

(n=10)


Sample 1 (25.2 U/ml): 4.2%

### Interassay-Precision:

(n=12)

Sample 1 (27.6 U/ml): 3.8%

## Typical Standard Curve



## Assay Procedure

Incubation of samples/ standards/ controls	100 µl	2 hrs, 4°C
Wash		
Incubation of detection antibody	100 µl	1 hr, room temperature
Wash		
Substrate incubation	100 µl/well	20 min, room temperature
Add Stopp solution	100 µl/well	
Read at 450nm		

## Order informations

Product	Catalog number	Price (€)
EIA for Quantitative Determination of anti-β2 AR AA, 1x96 determ.	12700	950.-

Ⓒ In vitro-Diagnostikum